

A18151

Leader in Biomolecular Solutions for Life Science



Cre Rabbit pAb

Catalog No.: A18151

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

ALL

Background

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-343 of e. coli Cre (WP_097745291.1)

Synonyms

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.