

A18201

Leader in Biomolecular Solutions for Life Science



pan-Ubiquitin Remnant Motif (K-ε-GG) Rabbit pAb

Catalog No.: A18201

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

ALL

Background

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

Swiss Prot

Immunogen

A synthetic peptide corresponding to a sequence containing ubiquitinated K.

Synonyms

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.