

A18269

Leader in Biomolecular Solutions for Life Science



NAD7 Rabbit pAb

Catalog No.: A18269

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Zea mays

Background

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

4055893

Swiss Prot

Immunogen

Recombinant protein of Zea mays NAD7.

Synonyms

NAD7

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.