

ACP6 Rabbit pAb

Catalog No.: A18456

Basic Information

Observed MW

Refer to figures

Calculated MW

49kDa

Category

Primary antibody

Applications

ELISA, IHC-P

Cross-Reactivity

Human, Mouse, Rat

Background

This gene encodes a member of the histidine acid phosphatase protein family. The encoded protein hydrolyzes lysophosphatidic acid, which is involved in G protein-coupled receptor signaling, lipid raft modulation, and in balancing lipid composition within the cell. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

IHC-P 1:50 - 1:200

Immunogen Information

Gene ID

51205

Swiss Prot

Q9NPH0

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-428 of human ACP6 (NP_057445.4).

Synonyms

LPAP; ACPL1; PACPL1; ACP6

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.