

ACOT4 Rabbit pAb

Catalog No.: A18535

Basic Information

Observed MW

Refer to figures

Calculated MW

46kDa

Category

Primary antibody

Applications

ELISA, IHC-P

Cross-Reactivity

Human, Mouse, Rat

Background

Enables acyl-CoA hydrolase activity and succinyl-CoA hydrolase activity. Involved in carboxylic acid metabolic process; saturated monocarboxylic acid metabolic process; and succinyl-CoA metabolic process. Located in peroxisome.

Recommended Dilutions

IHC-P 1:50 - 1:200

Immunogen Information

Gene ID

122970

Swiss Prot

Q8N9L9

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 180-340 of human ACOT4 (NP_689544.3).

Synonyms

PTE1B; PTE2B; PTE-Ib; ACOT4

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.