

A18679

Leader in Biomolecular Solutions for Life Science

LOC101736823 Rabbit pAb



Catalog No.: A18679 **2 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Bombyx mori

Background

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID
101736823

Swiss Prot

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 251-460 of Bombyx mori LOC101736823.

Synonyms

SQSTM1; LOC101736823

Contact



www.abclonal.com

Product Information

Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.