

# Methyl-eIF4A1-H132 Rabbit pAb

Catalog No.: A18814

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

46kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Human

## Background

---

Enables double-stranded RNA binding activity. Predicted to be involved in cytoplasmic translational initiation. Located in cytoplasm.

## Recommended Dilutions

---

WB 1:500 - 1:2000

## Immunogen Information

---

### Gene ID

1973

### Swiss Prot

P60842

### Immunogen

A synthetic peptide of human Methyl-eIF4A1-H132, 3-methyl (tau-methyl).

### Synonyms

DDX2A; EIF4A; EIF-4A; eIF4A-I; eIF-4A-I; Methyl-eIF4A1-H132

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.