

A18843

Leader in Biomolecular Solutions for Life Science



DNA-RNA Hybrid Mouse mAb

Catalog No.: A18843

Basic Information

Observed MW

Calculated MW

Category

Small Molecule-specific Antibody

Applications

ELISA,DB

Cross-Reactivity

Human,Mouse,Rat,Other (Wide Range Predicted)

CloneNo number

AMC0335

Background

Recommended Dilutions

DB 1:500 - 1:2000

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Synthetic DNA/RNA

Synonyms

Contact



www.abclonal.com

Product Information

Source

Mouse

Isotype

IgG2a k

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.