

A18871

Leader in Biomolecular Solutions for Life Science



# N6,2'-O-dimethyladenosine/m6Am Rabbit pAb

Catalog No.: A18871

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

### Category

Small Molecule-specific Antibody

### Applications

ELISA,DB

### Cross-Reactivity

Human,Mouse,Rat,Other (Wide Range Predicted)

## Background

---

## Recommended Dilutions

---

DB 1:500 - 1:2000

## Immunogen Information

---

Gene ID

Swiss Prot

### Immunogen

Chemical compounds corresponding to N6,2'-O-dimethyladenosine / m6Am.

### Synonyms

## Contact

---



[www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.