

pan-Phospho-(Ser/Thr) ATM/ATR Substrate Rabbit pAb

Catalog No.: A19337

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

WB

Cross-Reactivity

All

Background

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

无

Swiss Prot

无

Immunogen

A synthetic phosphorylated peptide around pan Ser/Thr of ATM/ATR Substrate .

Synonyms

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.