

A19511

Leader in Biomolecular Solutions for Life Science



AHI1 Rabbit pAb

Catalog No.: A19511

Basic Information

Observed MW

Refer to figures

Calculated MW

137kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human

Background

This gene is apparently required for both cerebellar and cortical development in humans. This gene mutations cause specific forms of Joubert syndrome-related disorders. Joubert syndrome (JS) is a recessively inherited developmental brain disorder with several identified causative chromosomal loci. Alternatively spliced transcript variants encoding different isoforms have been identified.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

54806

Swiss Prot

Q8N157

Immunogen

A synthetic Peptide of human AHI1

Synonyms

ORF1; AHI-1; JBTS3; dj71N10.1; AHI1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, pH7.3.