

Rabbit anti GFP-Tag mAb

Catalog No.: A19535 **3 Publications**

Basic Information

Observed MW

26kDa

Calculated MW

26kDa

Category

Primary antibody

Applications

ELISA, WB, IF/ICC

Cross-Reactivity

ALL

CloneNo number

ARC2165

Recommended Dilutions

WB 1:2000 - 1:6000

IF/ICC 1:50 - 1:200

Background

Immunogen Information

Gene ID

Swiss Prot

无

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-239 of GFP.

Synonyms

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of extracts of normal Eukaryotic expression of GFP and Eukaryotic expression of GFP transfected with GFP Protein, using GFP antibody (A19535) at 1:5000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 1s.