

A19588

Leader in Biomolecular Solutions for Life Science



# CTCF Rabbit mAb

Catalog No.: A19588 **Recombinant** **1 Publications**

## Basic Information

### Observed MW

140kDa

### Calculated MW

83kDa

### Category

SMab Recombinant Monoclonal Antibody

### Applications

WB,IF/ICC,IP,ChIP,ChIP-seq,ELISA,CUT&Tag

### Cross-Reactivity

Human,Mouse,Rat

### CloneNo number

ARC0067

## Background

This gene is a member of the BORIS + CTCF gene family and encodes a transcriptional regulator protein with 11 highly conserved zinc finger (ZF) domains. This nuclear protein is able to use different combinations of the ZF domains to bind different DNA target sequences and proteins. Depending upon the context of the site, the protein can bind a histone acetyltransferase (HAT)-containing complex and function as a transcriptional activator or bind a histone deacetylase (HDAC)-containing complex and function as a transcriptional repressor. If the protein is bound to a transcriptional insulator element, it can block communication between enhancers and upstream promoters, thereby regulating imprinted expression. Mutations in this gene have been associated with invasive breast cancers, prostate cancers, and Wilms' tumors. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

## Recommended Dilutions

**WB** 1:500 - 1:2000

**IF/ICC** 1:50 - 1:200

**IP** 0.5µg-4µg antibody for 200µg-400µg extracts of whole cells

**ChIP** 3µg antibody for 10µg-15µg of Chromatin

**CUT&Tag** 10<sup>5</sup> cells /2 µg

**ChIP-seq** 1:50 - 1:100

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Immunogen Information

### Gene ID

10664

### Swiss Prot

P49711

### Immunogen

A synthetic peptide corresponding to a sequence within amino acids 628-727 of human CTCF (P49711).

### Synonyms

MRD21; FAP108; CFAP108; CTCF

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification


Affinity purification

### Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.

## Validation Data





CUT&Tag was performed using the CUT&Tag Assay Kit (pAG-Tn5) for Illumina (RK20265) from  $10^5$  K562 cells with 2  $\mu$ g CTCF Rabbit mAb along with a Goat Anti-Rabbit IgG(H+L). The CUT&Tag results indicate the enrichment pattern of CTCF in representative gene loci (MYC), as shown in figure.



Chromatin immunoprecipitations were performed with cross-linked chromatin from 293T cells and CTCF Rabbit mAb (A19588). The ChIP sequencing results indicate the enrichment pattern of CTCF in selected genomic region and representative gene loci (MYC), as shown in figure.


Western blot analysis of various lysates using CTCF Rabbit mAb (A19588) at 1:1000 dilution.  
 Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
 Lysates/proteins: 25 $\mu$ g per lane.  
 Blocking buffer: 3% nonfat dry milk in TBST.  
 Detection: ECL Basic Kit (RM00020).  
 Exposure time: 30s.



Confocal imaging of U-2 OS cells using CTCF

## Validation Data

Rabbit mAb (A19588, dilution 1:200) followed by a further incubation with Cy3 Goat Anti-Rabbit IgG (H+L) (AS007, dilution 1:500) (Red). The cells were counterstained with  $\alpha$ -Tubulin Mouse mAb (AC012, dilution 1:400) followed by incubation with ABflo® 488-conjugated Goat Anti-Mouse IgG (H+L) Ab (AS076, dilution 1:500) (Green). DAPI was used for nuclear staining (Blue). Objective: 100x.



Immunoprecipitation analysis of 300  $\mu$ g extracts of 293T cells using 3  $\mu$ g CTCF antibody (A19588). Western blot was performed from the immunoprecipitate using CTCF antibody (A19588) at a dilution of 1:1000.

