

A20402

Leader in Biomolecular Solutions for Life Science



# PbrGEF11 Rabbit pAb

Catalog No.: A20402

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

### Category

Polyclonal Antibody

### Applications

WB, ELISA

### Cross-Reactivity

Pyrus bretschneideri

## Background

---

## Recommended Dilutions

---

WB 1:500 - 1:2000

## Immunogen Information

---

### Gene ID

103962981

### Swiss Prot

### Immunogen

A synthetic peptide of Pyrus bretschneideri PbrGEF11.

### Synonyms

## Contact

---



[www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.