

A20403

Leader in Biomolecular Solutions for Life Science



PbrPOE1 Rabbit pAb

Catalog No.: A20403

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Pyrus bretschneideri

Background

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID
103959040

Swiss Prot

Immunogen

A synthetic peptide of Pyrus bretschneideri PbrPOE1.

Synonyms

Contact



www.abclonal.com

Product Information

Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.