

A20404

Leader in Biomolecular Solutions for Life Science



Dicer1 Rabbit pAb

Catalog No.: A20404

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Gallus gallus domesticus

Background

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

423437

Swiss Prot

Immunogen

Recombinant protein of Gallus gallus domesticus Dicer1.

Synonyms

dicer; Dicer1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.