

ADAM20 Rabbit pAb

Catalog No.: A20442

Basic Information

Observed MW

Refer to figures

Calculated MW

82kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. The expression of this gene is testis-specific.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

8748

Swiss Prot

O43506

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 257-470 of human ADAM20 (NP_003805.3).

Synonyms

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.