

A20586

Leader in Biomolecular Solutions for Life Science



# OsWRKY67 Rabbit pAb

Catalog No.: A20586

## Basic Information

---

### Observed MW

25kDa

### Calculated MW

21kDa

### Category

Polyclonal Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Oryza sativa

## Background

---

WRKY genes encode transcription factors that are involved in the regulation of various biological processes

## Recommended Dilutions

---

WB 1:500 - 1:2000

## Immunogen Information

---

### Gene ID

4337998

### Swiss Prot

Q65WY5

### Immunogen

Recombinant protein of Oryza sativa OsWRKY67

### Synonyms

OsWRKY67

## Contact

---



[www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

## Validation Data

---



Western blot analysis of extracts of various tissues from the japonica rice (*Oryza sativa* L.) variety Zhonghua 11, using OsWRKY67 Rabbit pAb (A20586) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Enhanced Kit (RM00021).  
Exposure time: 120s.