

A20622

Leader in Biomolecular Solutions for Life Science



NOG1 Rabbit pAb

Catalog No.: A20622

Basic Information

Observed MW

43kDa

Calculated MW

43kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Oryza sativa

Background

encodes an enoyl-CoA hydratase/isomerase, increases the grain yield of rice by enhancing grain number per panicle without a negative effect on the number of panicles per plant or grain weight.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

4324637

Swiss Prot

Q5JMX9

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of Oryza sativa NOG1. (Q5JMX9).

Synonyms

NOG1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.

Validation Data



Western blot analysis of extracts of various tissues from the japonica rice (*Oryza sativa* L.) variety Zhonghua 11, using NOG1 Rabbit pAb (A20622) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Enhanced Kit (RM00021).
Exposure time: 300s.