

A20641

Leader in Biomolecular Solutions for Life Science



GL7 Rabbit pAb

Catalog No.: A20641

Basic Information

Observed MW

104kDa

Calculated MW

104kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Oryza sativa

Background

CTB4a (cold tolerance at booting stage), encoding a conserved leucine-rich repeat receptor-like kinase. The different CTB4a alleles confer distinct levels of cold tolerance and selection for variation in the CTB4a promoter region has occurred on the basis of environmental temperature.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

4343837

Swiss Prot

Q6YW00

Immunogen

Recombinant protein of Oryza sativa GL7.

Synonyms

GL7

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.

Validation Data



Western blot analysis of extracts of various tissues from the japonica rice (*Oryza sativa* L.) variety Zhonghua 11, using GL7 Rabbit pAb (A20641) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Enhanced Kit (RM00021).
Exposure time: 30s.