

A2425

Leader in Biomolecular Solutions for Life Science



## WWP2 Rabbit pAb

Catalog No.: A2425

1 Publications

### Basic Information

#### Observed MW

105kDa

#### Calculated MW

99kDa

#### Category

Polyclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Human

### Background

This gene encodes a member of the Nedd4 family of E3 ligases, which play an important role in protein ubiquitination. The encoded protein contains four WW domains and may play a role in multiple processes including chondrogenesis and the regulation of oncogenic signaling pathways via interactions with Smad proteins and the tumor suppressor PTEN. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 10.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

#### Gene ID

11060

#### Swiss Prot

O00308

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 151-400 of human WWP2 (NP\_001257383.1).

#### Synonyms

AIP2; WWp2-like; WWP2

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification


Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of lysates from A-549 cells, using WWP2 Rabbit pAb (A2425).