

A3389

Leader in Biomolecular Solutions for Life Science



# FAM120A Mouse mAb

Catalog No.: A3389

## Basic Information

---

### Observed MW

122kDa

### Calculated MW

122kDa

### Category

SMab Recombinant Monoclonal Antibody

### Applications

WB, IHC-P, IF/ICC, ELISA

### Cross-Reactivity

Human

## Background

---

Enables RNA binding activity. Located in cytosol.

## Recommended Dilutions

---

WB	1:500 - 1:1000
IHC-P	1:20 - 1:50
IF/ICC	1:20 - 1:50

## Immunogen Information

---

### Gene ID

23196

### Swiss Prot

Q9NZB2

### Immunogen

A synthetic peptide of human FAM120A

### Synonyms

OSSA; C9orf10; HBVPTPAP; FAM120A

## Contact

---



[www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Mouse

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.  
Buffer: PBS with 0.02% sodium azide, pH 7.3.