

KANK2 Rabbit pAb

Catalog No.: A3406

Basic Information

Observed MW

100kDa

Calculated MW

91kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human, Mouse

Background

This gene encodes a member of the KN motif and ankyrin repeat domains (KANK) family of proteins, which play a role in cytoskeletal formation by regulating actin polymerization. The encoded protein functions in the sequestration of steroid receptor coactivators and possibly other proteins. Mutations in this gene are associated with impaired kidney podocyte function and nephrotic syndrome, and keratoderma and woolly hair.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

25959

Swiss Prot

Q63ZY3

Immunogen

A synthetic peptide of human KANK2

Synonyms

SIP; MXRA3; PPKWH; NPHS16; ANKRD25; KANK2

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.