

ACTR5 Rabbit pAb

Catalog No.: A3505

Basic Information

Observed MW

75kDa

Calculated MW

68kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human, Mouse, Hamster

Background

Involved in UV-damage excision repair and double-strand break repair. Located in cytoplasm and nucleus. Part of Ino80 complex.

Recommended Dilutions

WB 1:1000 - 1:2000

Immunogen Information

Gene ID

79913

Swiss Prot

Q9H9F9

Immunogen

A synthetic peptide of human ACTR5

Synonyms

Arp5; INO80M; ACTR5

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using ACTR5 antibody (A3505).
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.