

# PCNXL2 Rabbit pAb

Catalog No.: A3510

## Basic Information

---

### Observed MW

253kDa

### Calculated MW

237kDa

### Category

Primary antibody

### Applications

WB

### Cross-Reactivity

Human

## Background

---

This gene contains coding mononucleotide repeats that are associated with tumors of high microsatellite instability (MSI-H). Defects in this gene are involved in the tumorigenesis of MSI-H colorectal carcinomas.

## Recommended Dilutions

---

WB 1:500 - 1:1000

## Immunogen Information

---

### Gene ID

80003

### Swiss Prot

A6NKB5

### Immunogen

A synthetic peptide of human PCNXL2

### Synonyms

PCNXL2

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.