

ZKSCAN3 Rabbit pAb

Catalog No.: A3514

Basic Information

Observed MW

61kDa

Calculated MW

61kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human

Background

Enables DNA-binding transcription repressor activity, RNA polymerase II-specific; RNA polymerase II cis-regulatory region sequence-specific DNA binding activity; and chromatin binding activity. Involved in several processes, including negative regulation of autophagy; negative regulation of cellular senescence; and regulation of transcription, DNA-templated. Located in cytoplasm and nucleoplasm.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

80317

Swiss Prot

Q9BRR0

Immunogen

A synthetic peptide of human ZKSCAN3

Synonyms

ZF47; Zfp47; ZFP306; ZNF306; ZNF309; zfp-47; ZSCAN13; ZSCAN35; dj874C20.1; dj874C20.1.; ZKSCAN3

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, pH7.3.