

A3532

Leader in Biomolecular Solutions for Life Science



ABHD14B Rabbit pAb

Catalog No.: A3532

Basic Information

Observed MW

22kDa

Calculated MW

22kDa

Category

Mouse Monoclonal Antibody

Applications

WB,FC

Cross-Reactivity

Human

Background

Enables hydrolase activity. Involved in positive regulation of transcription by RNA polymerase II. Located in cytosol; nucleolus; and nucleoplasm.

Recommended Dilutions

WB 1:500 - 1:2000

FC 1:20 - 1:50

Immunogen Information

Gene ID

84836

Swiss Prot

Q96IU4

Immunogen

A synthetic peptide of human ABHD14B

Synonyms

CIB; HEL-S-299; ABHD14B

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG


Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, pH7.3.

Validation Data



Western blot analysis of extracts of LO2 cells, using ABHD14B antibody (A3532).
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.