

GABPB2 Rabbit pAb

Catalog No.: A3543

Basic Information

Observed MW

55kDa

Calculated MW

49kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Mouse, Rat

Background

Enables transcription cis-regulatory region binding activity. Acts upstream of or within positive regulation of transcription by RNA polymerase II. Located in nucleus.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

126626

Swiss Prot

Q8TAK5

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 239-448 of human GABPB2 (NP_653219.1).

Synonyms

GABPB-2; GABPB2

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using GABPB2 antibody (A3543) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 30s.