

DCAF4L2 Rabbit pAb

Catalog No.: A3547

Basic Information

Observed MW

Refer to figures

Calculated MW

44kDa

Category

Primary antibody

Applications

WB,FC

Cross-Reactivity

Human, Mouse, Rat

Background

Recommended Dilutions

WB	1:500 - 1:2000
FC	1:20 - 1:50

Immunogen Information

Gene ID	Swiss Prot
138009	Q8NA75

Immunogen

A synthetic Peptide of human DCAF4L2

Synonyms

WDR21C; DCAF4L2

Contact

 | www.abclonal.com

Product Information

Source	Isotype	Purification
Rabbit	IgG	Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide,pH7.3.