

TMF1 Rabbit pAb

Catalog No.: A3642

Basic Information

Observed MW

131kDa

Calculated MW

123kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human

Background

Enables androgen receptor binding activity and nuclear receptor coactivator activity. Involved in androgen receptor signaling pathway and positive regulation of transcription by RNA polymerase II. Located in Golgi apparatus.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

7110

Swiss Prot

P82094

Immunogen

A synthetic peptide of human TMF1

Synonyms

TMF; ARA160; TMF1

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, pH7.3.