

TBCB Rabbit pAb

Catalog No.: A3789

Basic Information

Observed MW

36kDa

Calculated MW

27kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human

Background

Predicted to be involved in cell differentiation and nervous system development.
Located in cytosol.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

1155

Swiss Prot

Q99426

Immunogen

A synthetic peptide of human TBCB

Synonyms

CG22; CKAP1; CKAPI; TBCB

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.

Validation Data



Western blot analysis of extracts of HeLa cells, using TBCB antibody (A3789).
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25 μ g per lane.
Blocking buffer: 3% nonfat dry milk in TBST.