

GALK2 Rabbit pAb

Catalog No.: A3875

Basic Information

Observed MW

56kDa

Calculated MW

50kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human

Background

This gene encodes a highly efficient N-acetylgalactosamine (GalNAc) kinase, which has galactokinase activity when galactose is present at high concentrations. The encoded protein is a member of the GHMP kinase family. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

2585

Swiss Prot

Q01415

Immunogen

A synthetic peptide of human GALK2

Synonyms

GK2; GALK2

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.

Validation Data



Western blot analysis of extracts of NCI-H460 cells, using GALK2 antibody (A3875).
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.