

A4291

Leader in Biomolecular Solutions for Life Science



EFTUD2 Rabbit pAb

Catalog No.: A4291

Basic Information

Observed MW

110kDa

Calculated MW

109kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human

Background

This gene encodes a GTPase which is a component of the spliceosome complex which processes precursor mRNAs to produce mature mRNAs. Mutations in this gene are associated with mandibulofacial dysostosis with microcephaly. Multiple transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

9343

Swiss Prot

Q15029

Immunogen

A synthetic peptide of human EFTUD2

Synonyms

MFDM; MFDGA; Snu114; Snrp116; SNRNP116; U5-116KD; EFTUD2

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.

Validation Data



Western blot analysis of extracts of HeLa cells, using EFTUD2 antibody (A4291).
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25 μ g per lane.
Blocking buffer: 3% nonfat dry milk in TBST.