

STARD8 Rabbit pAb

Catalog No.: A4336

Basic Information

Observed MW

125kDa

Calculated MW

113kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human

Background

This gene encodes a member of a subfamily of Rho GTPase activating proteins that contain a steroidogenic acute regulatory protein related lipid transfer domain. The encoded protein localizes to focal adhesions and may be involved in regulating cell morphology. This protein may also function as a tumor suppressor.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

9754

Swiss Prot

Q92502

Immunogen

A synthetic Peptide of human STARD8

Synonyms

DLC3; ARHGAP38; STARTGAP3; STARD8

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.