

SRGAP3 Rabbit pAb

Catalog No.: A4354

Basic Information

Observed MW

121kDa

Calculated MW

125kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human

Background

Predicted to enable GTPase activator activity. Predicted to be involved in negative regulation of cell migration. Predicted to be located in cytosol. Predicted to be active in cytoplasm.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

9901

Swiss Prot

O43295

Immunogen

A synthetic Peptide of human SRGAP3

Synonyms

WRP; MEGAP; SRGAP2; ARHGAP14; SRGAP3

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.