

ANGPTL7 Rabbit pAb

Catalog No.: A4386

Basic Information

Observed MW

Refer to figures

Calculated MW

40kDa

Category

Primary antibody

Applications

WB,FC

Cross-Reactivity

Human, Mouse

Background

Enables identical protein binding activity. Involved in negative regulation of vasculature development involved in avascular cornea development in camera-type eye and regulation of extracellular matrix organization. Located in extracellular region.

Recommended Dilutions

WB	1:500 - 1:2000
FC	1:20 - 1:50

Immunogen Information

Gene ID	Swiss Prot
10218	043827

Immunogen

A synthetic Peptide of human ANGPTL7

Synonyms

AngX; CDT6; dj647M16.1; ANGPTL7

Contact

 | www.abclonal.com

Product Information

Source	Isotype	Purification
Rabbit	IgG	Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide,pH7.3.