

# TNIP1 Rabbit pAb

Catalog No.: A4404

## Basic Information

---

### Observed MW

70kDa

### Calculated MW

72kDa

### Category

Primary antibody

### Applications

WB

### Cross-Reactivity

Human

## Background

---

This gene encodes an A20-binding protein which plays a role in autoimmunity and tissue homeostasis through the regulation of nuclear factor kappa-B activation. Mutations in this gene have been associated with psoriatic arthritis, rheumatoid arthritis, and systemic lupus erythematosus. Multiple transcript variants encoding different isoforms have been found for this gene.

## Recommended Dilutions

---

WB 1:500 - 1:2000

## Immunogen Information

---

### Gene ID

10318

### Swiss Prot

Q15025

### Immunogen

A synthetic peptide of human TNIP1

### Synonyms

VAN; NAF1; ABIN-1; nip40-1; TNIP1

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.

## Validation Data

---



Western blot analysis of extracts of HeLa cells, using TNIP1 antibody (A4404).  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.