

A4408

Leader in Biomolecular Solutions for Life Science



# HMG20B Rabbit pAb

Catalog No.: A4408

## Basic Information

---

### Observed MW

36kDa

### Calculated MW

36kDa

### Category

Mouse Monoclonal Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Human,Mouse

## Background

---

Predicted to enable DNA binding activity. Predicted to be involved in regulation of gene expression. Predicted to act upstream of or within negative regulation of protein sumoylation; positive regulation of neuron differentiation; and skeletal muscle cell differentiation. Located in nuclear body.

## Recommended Dilutions

---

WB 1:500 - 1:1000

## Immunogen Information

---

### Gene ID

10362

### Swiss Prot

Q9P0W2

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-317 of human HMG20B (NP\_006330.2).

### Synonyms

SOXL; HMGX2; BRAF25; BRAF35; HMGXB2; PP7706; pp8857; SMARCE1r; HMG20B

## Contact

---



[www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

## Validation Data

---



Western blot analysis of various lysates using HMG20B Rabbit pAb (A4408) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates / proteins: 25 µg per lane.  
Blocking buffer: 3 % nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 5s.