

A4492

Leader in Biomolecular Solutions for Life Science



ANXA10 Rabbit pAb

Catalog No.: A4492

Basic Information

Observed MW

37kDa

Calculated MW

37kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. The function of this gene has not yet been determined.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

11199

Swiss Prot

Q9UJ72

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-324 of human ANXA10 (NP_009124.2).

Synonyms

ANX14; ANXA10

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using ANXA10 antibody (A4492) at 1:3000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 10s.