

A4512

Leader in Biomolecular Solutions for Life Science



SCAF8 Rabbit pAb

Catalog No.: A4512

Basic Information

Observed MW

141kDa

Calculated MW

141kDa

Category

Mouse Monoclonal Antibody

Applications

WB

Cross-Reactivity

Human

Background

Enables RNA binding activity and RNA polymerase II C-terminal domain phosphoserine binding activity. Involved in negative regulation of termination of RNA polymerase II transcription, poly(A)-coupled and positive regulation of DNA-templated transcription, elongation. Located in nucleoplasm.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

22828

Swiss Prot

Q9UPN6

Immunogen

A synthetic Peptide of human SCAF8

Synonyms

RBM16; SCAF8

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, pH7.3.