

PDS5A Rabbit pAb

Catalog No.: A4537

Basic Information

Observed MW

Refer to figures

Calculated MW

151kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Rat

Background

The protein encoded by this gene binds to the cohesin complex and associates with chromatin through most of the cell cycle. The encoded protein may play a role in regulating sister chromatid cohesion during mitosis. Two transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

23244

Swiss Prot

Q29RF7

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1020-1120 of human PDS5A (Q29RF7).

Synonyms

PIG54; SCC112; SCC-112; PDS5A

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.