

# AGTPBP1 Rabbit pAb

Catalog No.: A4541

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

138kDa

### Category

Primary antibody

### Applications

WB

### Cross-Reactivity

Human, Mouse, Rat

## Background

---

NNA1 is a zinc carboxypeptidase that contains nuclear localization signals and an ATP/GTP-binding motif that was initially cloned from regenerating spinal cord neurons of the mouse.

## Recommended Dilutions

---

WB 1:200 - 1:2000

## Immunogen Information

---

### Gene ID

23287

### Swiss Prot

Q9UPW5

### Immunogen

A synthetic Peptide of human AGTPBP1

### Synonyms

CCP1; NNA1; CONDCA; AGTPBP1

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.