

HSPBP1 Rabbit pAb

Catalog No.: A4574

Basic Information

Observed MW

Refer to figures

Calculated MW

39kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

Enables ubiquitin protein ligase binding activity. Involved in positive regulation of proteasomal ubiquitin-dependent protein catabolic process and positive regulation of protein ubiquitination. Predicted to be active in endoplasmic reticulum.

Recommended Dilutions

WB 1:200 - 1:2000

Immunogen Information

Gene ID

23640

Swiss Prot

Q9NZL4

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-359 of human HSPBP1 (NP_036399.3).

Synonyms

FES1; HSPBP1

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.