

CNNM3 Rabbit pAb

Catalog No.: A4625

Basic Information

Observed MW

Refer to figures

Calculated MW

76kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

Predicted to enable transmembrane transporter activity. Predicted to be involved in ion transport; magnesium ion homeostasis; and transmembrane transport. Located in membrane.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

26505

Swiss Prot

Q8NE01

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 410-640 of human CNNM3 (NP_951060.1).

Synonyms

ACDP3; CNNM3

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.