

# UBE2S Rabbit pAb

Catalog No.: A4658 **1 Publications**

## Basic Information

### Observed MW

24kDa

### Calculated MW

24kDa

### Category

Primary antibody

### Applications

ELISA, WB, IF/ICC

### Cross-Reactivity

Human, Mouse, Rat

## Background

This gene encodes a member of the ubiquitin-conjugating enzyme family. The encoded protein is able to form a thiol ester linkage with ubiquitin in a ubiquitin activating enzyme-dependent manner, a characteristic property of ubiquitin carrier proteins.

## Recommended Dilutions

|               |                |
|---------------|----------------|
| <b>WB</b>     | 1:200 - 1:2000 |
| <b>IF/ICC</b> | 1:50 - 1:200   |

## Immunogen Information

|                |                   |
|----------------|-------------------|
| <b>Gene ID</b> | <b>Swiss Prot</b> |
| 27338          | Q16763            |

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-222 of human UBE2S (NP\_055316.2).

### Synonyms

EPF5; E2EPF; E2-EPF; UBE2S

## Contact

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

|               |                |                       |
|---------------|----------------|-----------------------|
| <b>Source</b> | <b>Isotype</b> | <b>Purification</b>   |
| Rabbit        | IgG            | Affinity purification |

### Storage

Store at -20°C. Avoid freeze / thaw cycles.  
Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.