

SPCS1 Rabbit pAb

Catalog No.: A4661

Basic Information

Observed MW

13kDa

Calculated MW

18kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human

Background

Predicted to enable peptidase activity and ribosome binding activity. Involved in viral protein processing and virion assembly. Is integral component of endoplasmic reticulum membrane.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

28972

Swiss Prot

Q9Y6A9

Immunogen

A synthetic Peptide of human SPCS1

Synonyms

SPC1; SPC12; HSPC033; YJR010C-A; SPCS1

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.