

A4786

Leader in Biomolecular Solutions for Life Science



RHOF Rabbit pAb

Catalog No.: A4786

Basic Information

Observed MW

27kDa

Calculated MW

24kDa

Category

Mouse Monoclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse

Background

Predicted to enable GTP binding activity; GTPase activity; and protein kinase binding activity. Involved in actin filament organization. Located in extracellular exosome.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

54509

Swiss Prot

Q9HBH0

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 40-120 of human RHOF (NP_061907.2).

Synonyms

RIF; ARHF; RHOF

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

Validation Data



Western blot analysis of various lysates using RHO Rabbit pAb (A4786) at 1:3000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Enhanced Kit (RM00021).
Exposure time: 90s.