

# ATL2 Rabbit pAb

Catalog No.: A4925

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

66kDa

### Category

Primary antibody

### Applications

WB,FC

### Cross-Reactivity

Human, Mouse, Rat

## Background

---

Enables identical protein binding activity. Involved in Golgi organization; endoplasmic reticulum tubular network membrane organization; and protein homooligomerization. Located in endoplasmic reticulum tubular network membrane. Is integral component of membrane.

## Recommended Dilutions

---

<b>WB</b>	1:500 - 1:2000
<b>FC</b>	1:20 - 1:50

## Immunogen Information

---

<b>Gene ID</b>	<b>Swiss Prot</b>
64225	Q8NHH9

### Immunogen

A synthetic Peptide of human ATL2

### Synonyms

aip-2; ARL3IP2; ARL6IP2; atlastin2; ATL2

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

<b>Source</b>	<b>Isotype</b>	<b>Purification</b>
Rabbit	IgG	Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.  
Buffer: PBS with 0.02% sodium azide,pH7.3.