

UBE2Z Rabbit pAb

Catalog No.: A4951

Basic Information

Observed MW

38kDa

Calculated MW

38kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human

Background

This gene encodes an enzyme which ubiquitinates proteins which participate in signaling pathways and apoptosis.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

65264

Swiss Prot

Q9H832

Immunogen

A synthetic Peptide of human UBE2Z

Synonyms

USE1; HOYS7; UBE2Z

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.

Validation Data



Western blot analysis of extracts of MDA-MB-453 cells, using UBE2Z antibody (A4951).
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.