

FAM118B Rabbit pAb

Catalog No.: A4970

Basic Information

Observed MW

54kDa

Calculated MW

39kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human

Background

Enables identical protein binding activity. Predicted to be involved in Cajal body organization. Predicted to be located in Cajal body.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

79607

Swiss Prot

Q9BPY3

Immunogen

A synthetic Peptide of human FAM118B

Synonyms

FAM118B

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.

Validation Data



Western blot analysis of extracts of Jurkat cells, using FAM118B antibody (A4970).
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.